



Reagent set, lead, dithizone method, (no air shipments)

Številka izdelka:

2243100

Nearno

ZASTAREL IZDELEK
Ta izdelek ni več na voljo.



For the determination of Lead, Dithizone Method.

For determination of lead by the Dithizone extraction method. USEPA accepted for reporting for wastewater analysis using Hach Method 8033. Range: 3 to 300 µg/L Pb. Set Includes reagents for 100 tests. For air freight shipments, order 2243101

Specifikacije

EPA compliant:	Da
Ime metode:	Dithizone
Inštrument:	All except DR/800, PCII
Merilno območje:	3 - 300 µg/L Pb
Shipping Method:	Lead Reagent set: No Air Freight
Številka metode:	8033
Število testov:	60 - 100