



Reagent set, aluminum, 50 mL sample, ECR method

Številka izdelka:

2443400

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nearno



For the determination of Aluminium, Eriochrome Cyanine R Method. For 25ml sample.

For determination of Aluminum by the Eriochrome Cyanine R (ECR) Method. Hach Method 8326. Sample Size: 50 mL Range: 0.002 - 0.220 mg/L. Reagent set contains ECR Reagent Powder Pillows, Hexamethylenetetramine Buffer Reagent Powder Pillows and ECR Masking Reagent Solution for approximately 100 tests.

Specifikacije

Ime metode:	Eriochrome Cyanine R
Inštrument:	DR/3000, DR/2000
Merilno območje:	0 - 0,220 mg/L Al
Številka metode:	8326
Število testov:	100