



Citrate Buffer Powder Pillows

Številka izdelka:

2107669

EUR Cena:

Kontakt

Na voljo

Nearno



Powder reagent in sealed packages for the determination of Manganese, cold periodate method. For 10 mL sample.

For determination of manganese by the periodate oxidation method. EPA approved for reporting wastewater analyses (digestion required). Hach method 8034. Sample size: 10 mL Pack of 100 Powder Pillows.

Specifikacije

Ime metode:	Periodate Oxidation
Platforma:	Prašasti reagenti
Prostornina / velikost paketa:	100 /pk
Število testov:	100