



Hardness (Calcium) Reagent Set, HR, Digital Titrator

Številka izdelka: 2447500
EUR Cena: Kontakt
Datum pošiljanja ni prikazan

Nearno



For the determination of Calcium Hardness, 100-4000 mg/L CaCO₃, EDTA Method.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikacije

Ime metode:	EDTA Titration
Merilno območje:	100 - 4000 mg/L CaCO ₃
Parameter:	Hardness (Calcium)
Platforma:	Digitalni Titrator
Posebna značilnost:	Reagent set, HR, Digital Titrator
Številka metode:	8204
Število testov:	40-100