



Detergents Reagent Set, 0.002-0.275 mg/L LAS (Air Freight Configuration)

Številka izdelka:

2446801

Nearno

ZASTAREL IZDELEK
Ta izdelek ni več na voljo.



For determination of Detergents (Surfactants) by the Crystal Violet method (Hach Method 8028). Range: 0.002 - 0.275 mg/L LAS. Reagent Set includes: Benzene, Buffer Solution and Detergent Reagent Powder Pillows for 75 tests. Note: Reagents are packaged for air freight shipment.

Specifikacije

Ime metode:	Crystal Violet
Merilno območje:	0.002 - 0.275 mg/L LAS
Število testov:	75