



## CDTA Digital titrator cartridge, 0.800 M

Številka izdelka:

1440301

Nearno

EUR Cena:

Kontakt

Na voljo



Cyclohexanediamine tetraacetic acid standard solution. For Calcium and Total Hardness determination by Digital Titrator titration. 0.800 M.  
Digits = mg/L CaCO<sub>3</sub> in 100 mL. Digits x 0.001 = mmole.

### Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

### Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

### Versatile

More than 50 parameter/range combinations are available

### Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

---

## Specifikacije

|                  |                                   |
|------------------|-----------------------------------|
| Ime metode:      | CDTA Titration                    |
| Merilno območje: | 100 - 4000 mg/L CaCO <sub>3</sub> |
| Platforma:       | Digitalni Titrator                |

---

## Obvezni pripomočki

- Zaščitne rokavice velikosti L, modre, nitril, brez pudra, 50 kosov (Item SM743L)