



Be Right™



Digital titrator cartridge, phenylarsine oxide (PAO), 0.00564 N

Številka izdelka: 199901
EUR Cena: Kontakt
Odprema v 2 tednih

Nearno



For low level chlorine determination by amperometric titration, using the Digital Titrator

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikacije

EPA compliant:	Da
Ime metode:	Amperometric, Forward
Merilno območje:	0 - 1000 µg/L Cl ₂
Parameter:	Chlorine, Free
Platforma:	Digitalni Titrator
Številka metode:	10026
Titant:	0.00564 N PAO