



Phenylarsine Oxide (PAO) Standard Solution, 0.00564 N, 1 L

Številka izdelka: 199953
EUR Cena:
Datum pošiljanja ni prikazan

Kontakt

Nearno



For chlorine determination by amperometric titration methods. Note: also available in Digital Titrator cartridges. 0.00564 N, 1 liter bottle

Specifikacije

Ime metode:	Amperometric Titration
Inštrument:	AutoCat [™] 9000
Merilno območje:	0,0012 - 5,000 mg/L Cl ₂
Platforma:	Raztopina
Številka metode:	8168