



EDTA digital titrator cartridge, 0.0716 M

Številka izdelka:

2081701

Nearno

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan



EDTA standard solution. For the determination of iron by Digital Titrator titration. 0.0716 M. Digits x 0.1 = mg/L Fe in 50 mL

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikacije

Ime metode:	EDTA
Merilno območje:	10 - 100 mg/L Fe
Parameter:	Iron
Platforma:	Digitalni Titrator
Število testov:	40 - 100
Titrant:	0.0716 M EDTA