



## MolyVer® 6 molybdenum reagent AccuVac® ampules, 25/pk

Številka izdelka:

2522025

Nearno

**ZASTAREL IZDELEK**  
Ta izdelek ni več na voljo.

Predlagani nadomestki:

- 2604100



For determination of HR Molybdate Molybdenum by the Mercaptoacetic Acid Method. Hach method 8036. Range: 0.2 - 40.0 mg/L Mo. Pack of 25 ampules. Also requires CDTA Solution, 0.4 M, Cat. No. 26154-36.

### Fast

AccuVac® Ampules provide quick results. Just snap, fill and read - saving you valuable time.

### Easy

No need to measure sample or reagents. The evacuated ampules draw in the correct sample volume, so anyone can obtain reliable results with minimal effort.

### Accurate

Optical quality evacuated vials assure precise measurements and reagent stability.

### Use them anywhere you need results.

Ampules are both reagent containers and photometric measurement cells. No need to transport reagents or samples.

---

## Specifikacije

|                  |                     |
|------------------|---------------------|
| Ime metode:      | Mercaptoacetic Acid |
| Merilno območje: | 0.2 - 40.0 mg/L Mo  |
| Parameter:       | Molybdenum          |
| Platforma:       | Accuvac             |
| Število testov:  | 25                  |