



Bottle, m-ColiBlue24 broth, 100 mL (50 test)

Številka izdelka: 2608442

EUR Cena: Kontakt

Odprema v 2 tednih

In 24 hours, Hach's m-ColiBlue24 Broth give you simultaneous results for total coliforms and E. coli. Hach's ready-to-use media eliminates measuring, mixing, and autoclaving steps necessary to prepare dehydrated media. Glass ampules (26084-20), plastic ampules (260845-00), 100mL bottles (26084-42) of m-ColiBlue24 and agar plates (28052-15), afford maximum shelf life and ease of use for pre-measured and prepared media. Enumerate total coliforms and E. coli on one petri dish. Differentiate colonies easily - red and blue indicate total coliforms and blue specifies E. coli. Get optimal recovery of stressed and injured organisms.

Use m-ColiBlue24 Broth to accurately monitor and evaluate total coliform and E. coli in drinking water; wastewater; bottled water; beverages; surface, ground, and well water; and ultrapure, chemical processing and pharmaceutical processing waters.

m-ColiBlue24 should only be used with pads manufactured by the following companies/brands: Pall Corporation (e.g., Sterile Petri Dishes with Absorbant Pads), Private Label/OEM products manufactured by Pall Gelman, Sartorius (cellulose and glass fiber pads), Sartorius BioSart products with the integrated Petri dish, and Millipore EO sterilized pads. No other pad is compatible with the product, and the use of any pad other than those recommended may result in unreliable, inaccurate data.

With ready-to-use m-ColiBlue24 Broth, you can:

Read and confirm results in 24 hours

Get superior sensitivity - 1 CFU/100mL

Conduct emergency testing and routine monitoring

Eliminate special equipment and ultraviolet lamp

Minimize background growth of non-coliform bacteria

Specifikacije

| | |
|--------------------------------|--------------------|
| Material: | Glass |
| Media Type: | m-ColiBlue24 Broth |
| Prostornina / velikost paketa: | 1 each |