



Reagent set, phosphorus, ascorbic acid rapid liquid™

Številka izdelka:

2678600

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nearno



For determination of reactive Phosphorus (ortho-phosphate) by the Ascorbic Acid Rapid Liquid™ method. Hach Method 10055. Range: 19 - 3000 µg/L PO₄. Set includes: Ascorbic Acid Reagent Dilution Solution, Ascorbic Acid Reagent Powder and Molybdate Reagent Solution for approximately 250 tests.

Convenience

Rapid Liquid methods utilize liquid reagents, bottletop dispensers and pour-thru cells to provide simplified testing.

Performance

Since the pour-thru cell provides a stationary optical system, Rapid Liquid methods typically provide better precision than methods that require movement of sample cells.

Efficiency

Rapid Liquid methods are ideal for high-volume testing.

Specifikacije

Ime metode:	Ascorbic Acid
Inštrument:	All except DR 900, DR/800 Series, PC II
Merilno območje:	19 - 3000 µg/L PO ₄
Platforma:	Raztopina
Številka metode:	10055
Število testov:	250