



## Free and total chlorine test strips, 0 - 10 mg/L, 50 tests

Številka izdelka: 2745050

EUR Cena: Kontakt

Odprema v 2 tednih

The most widely used disinfectant for drinking water, chlorine is also important for sanitizing swimming pools, cooling towers, other industrial equipment, and in the treatment of municipal wastewater. Its measurement and control are vital for both safety and economic reasons. Chlorine exists in water as "free" chlorine (hypochlorite) or "combined" chlorine (chloramines). These chlorine test strips streamline water testing for field technicians, enabling swift analysis of multiple samples within minutes and facilitating on-site evaluations without the need for extensive preparation or cleanup. Widely utilized in laboratories globally, our chlorine test strips serve as invaluable tools for pre-test screening, swiftly detecting potential interferences with lab testing procedures. Each package contains 50 tests.

---

### Specifikacije

Enota:	50 tests
Komentar:	*mg/L unless otherwise noted; ppb = $\mu\text{g/L}$ ; ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Merilno območje:	0 - 10 mg/L $\text{Cl}_2$
Najmanjši korak rasti:	Steps: 0, 0.5, 1, 2, 4, 10
Parameter:	Free and Total Chlorine
Platforma:	Test Strip
Test Kit model:	TS