



Flow through unit for NT3100 sc/NT3200 sc 5 mm, Nitratax plus sc 5 mm, Uvas plus sc 5 mm probe

Številka izdelka: LZX867

EUR Cena: Kontakt

Odprema v 1 tednu

Specifikacije

Delovna temperatura:	2 - 40 °C (sample temperature)
Dimenzija (v x š x g):	210 mm x 500 mm x 160 mm
Dolžina:	333 mm
Dolžina kabla:	10 m
Dolžina poti:	2 mm
Material:	Measuring cell: PVC;
	Glands: PVDF;
	Sample tube: PVC;
	Screws: Stainless Steel;
	Seals: EPDM
Materiali v stiku z vodo:	Measuring cell: PVC;
	Seals: EPDM;
	Glands: PVDF;
	Sample tube: PVC
Območje tlaka:	Not pressurized, open drain
Pretok vzorca:	Min. 0.5 L/h
Prostornina vzorca:	0,5 - 10 L
Temperatura vzorca:	2 - 40°C (35,6 - 104°F)
Teža:	3,6 kg
Vhod vzorca:	Tube ID 4 mm/AD 6 mm