



Sedimeter for NT3x00sc, Nitratax sc and Uvas sc

Številka izdelka: **LZX450**

EUR Cena: Kontakt

Odprema v 1 tednu

The sedimenter can be used as a bypass solution. The sample can be let through the sedimenter without prefiltering.

Specifikacije

Dolžina:	Intallation plate for wall mount: 800 mm x 495 mm
Inštrument:	Suitable for all NT3x00sc, Nitratax and Uvas sc sensors
Jeziki v navodilih za uporabo:	Non verbal instruction sheet
Material:	PVC
Območje tlaka:	Not pressurized, open drain
Teža:	Sedimeter incl. wall plate only: 14.7 kg