



Water, deionised and demineralised, 500mL

Številka izdelka:	27249
EUR Cena:	Kontakt
Na voljo	

Demineralised Type 2 Water

Specific Conductance: $\leq 1 \mu\text{S/cm}$,

Silica (SiO_2): $\leq 30 \mu\text{g/L}$,

Chloride (Cl^-): $\leq 20 \mu\text{g/L}$,

Hardness: Negative,

Consumption of KMnO_4 : Negative

Specifikacije

Obseg:	500 mL
Opis:	Water, Deionized
Platforma:	Raztopina