



**Be Right™**



## **Filling Solution, Internal, Polarographic DO, 30 mL**

**Številka izdelka:** 2759123

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Filling solution for Polarographic Dissolved Oxygen probe, Potassium Chloride/Glycerine, 30 mL bottle

---

### **Specifikacije**

|                                |                 |
|--------------------------------|-----------------|
| Formula:                       | KCl & Glycerine |
| Opis:                          | 5130, 5197000   |
| Prostornina / velikost paketa: | 30 mL           |
| Solution Type:                 | Refill          |