



Singlet Single-Use pH Buffer Kit, pH 4.01 & 7.00, pk/2x10

Številka izdelka: 2769920

EUR Cena: Kontakt

Odprema v 2 tednih

Singlet Solution Packs

pH Buffer Standard Solution for pH electrode calibration. Sealed foil pouches of pH buffers. pH is ± 0.02 at 25°C. Pack of 10 of pH 4.01 and 10 of pH 7.00, 20 mL each

Ideal for on-site analyses

Sealed pouches are easy to transport - no risk of spillage.

Easy to Use

Just open the pouch and insert the electrode.

Fresh and reliable

Foil pouches prevent contamination and air exposure.

Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

Specifikacije

Barva kodirana?:	Da
Koncentracija:	4,01 in 7,00 pH
Opis:	pH Buffer Calibration reagent
Package settings:	Singlet
pH vrednost(i):	4.01 & 7.00 ± 0.02 pH at 25°C
Pogoji skladiščenja:	10 - 25 °C (protect from light)
Posebni zapiski:	Discard solution after use
Prostornina / velikost paketa:	2 pks of 10 pcs
Sledljivo po NIST-u:	Da