



Standardna raztopina za prevodnost 1413 µS/cm; COA (analizni certifikat) je na voljo prek spleta

Številka izdelka: LZW9710.99
EUR Cena: Kontakt
Na voljo

Conductivity standard solutions for laboratory and field analyses

Standardna raztopina za prevodnost 1413 µS/cm; COA (analizni certifikat) je na voljo prek spleta o analizi; 250-mililitrska steklenička.

Wide selection of concentrations

and package sizes

Hach conductivity standards are offered in a large number of sizing to match each user environments and standard operating practices.

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard solutions that are used are key.

Design for use with Hach meters portfolio

The basis for accurate conductivity measurements is good care and calibration of the conductivity cell.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifikacije

Formula:	KCl
Koncentracija:	0.01M KCl
Opis:	Conductivity calibration solution
Package settings:	Bottle
Parameter:	Conductivity
Platforma:	Raztopina
Pogoji skladiščenja:	10 - 25 °C (protect from light)
Posebni zapiski:	Observed shelf life after opening, following recommended storage and use: 3 months
Prostornina / velikost paketa:	250 mL
Sledljivo po NIST-u:	Da