



**Be Right™**



## Test strips pH , 4-9 pH units

Številka izdelka: 2745650

EUR Cena: Kontakt

Na voljo

pH test strips offer a convenient and versatile solution for measuring the acidity or alkalinity of various aqueous solutions. Most natural waters exhibit pH levels ranging from pH 4 to pH 9, with common values exceeding pH 7 due to the presence of carbonates and bicarbonates, contributing to alkalinity. In industrial settings, maintaining precise pH levels in processes such as manufacturing, chemical production, and wastewater treatment is crucial for optimal efficiency and product quality. Similarly, in boiler feed water systems, controlling pH within a narrow range is essential to prevent corrosion and scale buildup. Additionally, pH regulation in swimming pools is vital for ensuring water comfort, swimmer safety, and the effectiveness of sanitizers.

With 50 tests per package, our pH test strips provide a cost-effective and user-friendly solution for accurate pH measurement in a wide range of applications.

**Easy to use, disposable, and inexpensive**

**Great for quick evaluations of pH levels in the lab or in the field**

---

### Specifikacije

Enota:	50 tests
Merilno območje:	4 - 9 pH units
Najmanjši korak rasti:	Steps: 4, 5, 6, 7, 8, 9 pH
Parameter:	pH, mid range
Platforma:	Test Strip
Test Kit model:	Test Strips