



Testkit, glycol visual, model EG-1

Številka izdelka:

2186400

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nearno



Semiquantitative determination of glycol in water or oil. Visual appearance of colour indicates glycol presence.

Ethylene glycol (anti-freeze) contamination of automotive lubricants is an indication of a cracked engine block. Glycol presence also can indicate leaks in cooling systems. Test Kit Model EG-1 presents a purple color if ethylene glycol and other 1,2 glycols are present in oil or water. The kit contains approximately 25 tests and comes with a 14.4" x 8" x 7" Blue Polypropylene Carrying Case.

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Ime metode:	Visual/
	Appearance of color
Merilno območje:	Presence or absence of glycol, 0- greater than 150 mg/L
Parameter:	Glycol - Ethylene glycol
	in oil or water
Število testov:	25
Test Kit model:	EG-1
Vrsta primera:	G