



TEST KIT, 5-EP mg/L HARDNESS

Številka izdelka:

145401

Nearno

EUR Cena:

Kontakt

Na voljo



Total Hardness drop count titration/EDTA test kit with range of 20 - 400 mg/L; contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Ime metode:	Drop count titration/ EDTA-powder pillows
Merilno območje:	20 - 400 mg/L CaCO ₃
Parameter:	Hardness, total - As CaCO ₃
Platforma:	Drop Count Titration
Smallest increment:	20
Število testov:	100
Test Kit model:	5-EP MG-L
Vrsta primera:	B