



Be Right™



## Hardness (Total & Calcium) Test Kit, Model HA-4P MG-L

Številka izdelka:

145701

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nevorno



Hardness, Total & Calcium (Magnesium by difference) drop count titration/EDTA test kit, range 20 - 400 mg/L CaCO<sub>3</sub>, 100 tests.

The Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field.

**Utilises drop count titration - a simple, economical titration method**

### Specifikacije

Ime metode:	Drop count titration/EDTA
Merilno območje:	20 - 400 mg/L CaCO <sub>3</sub>
Parameter:	Hardness, total & Calcium - (Magnesium by difference)
Platforma:	Drop Count Titration
Smallest increment:	20
Število testov:	100
Test Kit model:	HA-4P
	MG-L
Vrsta primera:	D
Vsebina paketa:	Includes five reagents, bottle, measuring tube, instruction sheet, and carrying case.

### Vsebina paketa

Includes five reagents, bottle, measuring tube, instruction sheet, and carrying case.