



Hardness (Total & Calcium) Test Kit, Model HAC-DT

Številka izdelka:

2063900

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nearno



Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10 - 4000 mg/L. Contains approximately 100 of each (Total and Calcium) test.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Specifikacije

Ime metode:	Digital Titrator/EDTA
Merilno območje:	10 - 4,000 mg/L CaCO_3
Parameter:	Hardness, total & Calcium - As CaCO_3
Platforma:	Digitalni Titrator
Smallest increment:	0.1-10
Število testov:	100 Ca
	+100 T
Test Kit model:	HAC-DT
Vrsta primera:	F
Vsebina paketa:	Includes 6 reagents, digital titrator, delivery tubes, graduated cylinder, flask, manual and carrying case.

Vsebina paketa

Includes 6 reagents, digital titrator, delivery tubes, graduated cylinder, flask, manual and carrying case.