



## Test kit, total chlorine, model CN-DT, 20-2000 mg/L

Številka izdelka: 2471100

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

**Nearno**



The most widely used disinfectant for drinking water, chlorine is also important for sanitizing swimming pools, cooling towers, other industrial equipment, and in the treatment of municipal wastewater. Its measurement and control are vital for both safety and economic reasons. Kit utilizes the Hach Digital Titrator and contains 100 tests.

**Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration**

**Rugged titration device is ideal for field use, with no glass burets required**

**Kit contains all required reagents and apparatus to perform analysis in the field**

### Specifikacije

Ime metode:	Digital Titrator/Thiosulfate
Komentar:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Merilno območje:	20 - 2,000 mg/L Cl <sub>2</sub>
Parameter:	Chlorine, mid and high range - As Total Cl <sub>2</sub>
Platforma:	Digitalni Titrator
Smallest increment:	0.2 - 5.0
Število testov:	100
Test Kit model:	CN-DT
Vrsta primera:	F