



Dissolved oxygen test kit, model OX-2P

Številka izdelka:

146900

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nearno



Low dissolved oxygen levels in surface water systems usually indicate serious pollution. Adequate amounts are crucial for fish life, but conversely, DO must be excluded from boiler feed water to prevent corrosion. This test kit utilizes drop count titration and contains 100 tests.

Method / Chemistry: Drop count titration

Replacement reagent set PN 2439802

Utilizes drop count titration - a simple, economical titration method

Specifikacije

Ime metode:	Drop count titration/ Modified Winkler
Merilno območje:	0.2 - 4.0 mg/L & 1 - 20 mg/L O ₂
Parameter:	Dissolved Oxygen - As DO
Platforma:	Drop Count Titration
Smallest increment:	0.2
Število testov:	100
Test Kit model:	OX-2P
Vrsta primera:	D