



Be Right™



Sulfite Test Kit, Model SU-5

Številka izdelka:

148002

EUR Cena:

Kontakt

Na voljo

Nevarno



Utilizing drop count titration/iodometric to determine sulfite, this kit contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Ime metode: Drop count,

titration Iodometric

Merilno območje: 1.0 - 20 mg/L & 10 - 200 mg/L SO₃²⁻

Parameter: Sulfite

Platforma: Drop Count Titration

Smallest increment: 1

10

Število testov: 100

Test Kit model: SU-5

Vrsta primera: D