



Test kit, sulfite, model SU-DT

Številka izdelka:

2063300

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nearno



Using the Digital Titrator, this kit contains 100 tests with a range of 10-800 mg/L.

Digital Titrator and standardized titration cartridges provide laboratory accuracy for field titrations

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Ime metode:	Digital Titrator/Iodometric
Merilno območje:	4 - 800 mg/L SO_3^{2-}
Najmanjši korak rasti:	0.4-4.0**
Parameter:	Sulfite - As SO_3^{2-}
Platforma:	Digitalni Titrator
Število testov:	100
Test Kit model:	SU-DT
Vrsta primera:	F