



Test kit, alkalinity, model AL-DT

Številka izdelka:

2063700

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nearno



Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Determine phenolphthalein and total alkalinity by visual endpoint titration. Using the Digital Titrator, this contains 100 tests with a range of 10-4,000 mg/L.

Digital Titrator and standardized titration cartridges provide laboratory accuracy for field titrations

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Ime metode:	Digital Titrator/Sulfuric acid
Komentar:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Merilno območje:	10 - 4,000 mg/L CaCO ₃
Najmanjši korak rasti:	0.1-10
Parameter:	Alkalinity - (P) & Total (MO) as CaCO ₃
Platforma:	Digitalni Titrator
Število testov:	100
Test Kit model:	AL-DT
Vrsta primera:	F