



**Be Right™**



## Test kit, alkalinity, model AL-TA

Številka izdelka:

2314500

Nevarno

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan



Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Utilizing drop count titration/sulfuric acid to determine alkalinity, this kit contains 100 tests.

**Utilizes drop count titration - a simple, economical titration method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

### Specifikacije

Ime metode:	Drop count titration/Sulfuric acid
Komentar:	*mg/L unless otherwise noted, ppb = $\mu\text{g}/\text{L}$ , ppm = mg/L; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Merilno območje:	23 - 495 gpg & 385 - 8500 mg/L
Parameter:	Alkalinity - (P) & Total (MO) As $\text{CaCO}_3$
Platforma:	Drop Count Titration
Smallest increment:	23 gpg
	385
Število testov:	100
Test Kit model:	AL-TA
Vrsta primera:	B