



Test kit, alkalinity, model AL-AP

Številka izdelka:

2444300

EUR Cena:

Kontakt

Odprema v 2 tednih

Nearno



Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Utilizing drop count titration/sulfuric acid to determine alkalinity, this kit contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Ime metode:	Drop count titration/Sulfuric acid
Komentar:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. **1 gpg = 17.1 mg/L (MO) = Methyl Orange
Merilno območje:	0 - 8 gpg & 1 - 20 gpg
Parameter:	Alkalinity - (P) & Total (MO) As CaCO ₃
Platforma:	Drop Count Titration
Smallest increment:	0.4 gpg
	1 gpg
Število testov:	100
Test Kit model:	AL-AP
Vrsta primera:	B