



BAG, WHIRL-PAK* 5-1/2X9 PK/10

Številka izdelka: 1437297

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Complies with Standard Methods and USEPA for drinking water sample collection. Made from poly resins approved by USDA and FDA for food contact. Sterile (certified), disposable, transparent polyethylene bags with trusted "Whirl-Pak" closure. Suitable for liquid or solid sample collection. "Write-on" surface accepts waterproof ink. 24 oz bag perfect for incubating petri dishes in water bath applications.

Specifikacije

Kapaciteta: 720 mL

Prostornina / velikost paketa: 10