



Be Right™



Paddle Tester, disinfection control

Številka izdelka:

2619510

EUR Cena:

Kontakt

Na voljo

Total Aerobic Bacteria/Disinfection Control

TSA/DE Neutralizing Agar, pk/10

Hach Paddle Testers allow fast, easy screening of water samples, solid surfaces, and nonviscous liquids where unsanitary conditions may exist. The paddle is a double-sided slide attached to the vial cap. Each side of the slide is used to perform a separate test. The slide has a molded-in grid for easy colony counting without opening the optically clear, leakproof vial.

Paddle Tester, Total Aerobic Bacteria/Disinfection Control (Tryptic Soy Agar/DE Neutralizing Agar), pk/10.

Specifikacije

Media Type: Total Aerobic Bacteria/Disinfection Control

Prostornina / velikost paketa: 10 /pk