



ToxTrak Toxicity Reagent Set, 25-49 tests

Številka izdelka:

2597200

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nearno



ToxTrak Toxicity Test Kit

Use the ToxTrak Toxicity Reagent Set to analyze your water or wastewater samples for toxicity. Follows Hach Toxicity Method 10017.

The ToxTrak reagents perform toxicity screens on drinking water, wastewater and natural water.

Toxicity is compared to a baseline using indigenous biomass or stock bacteria cultures with a range of 0 -100% inhibition.

A colorimetric method is used to determine the toxicity of samples, with analysis done using one of the Hach's trusted spectrophotometers.

Specifikacije

Inštrument:	All except DR/3000, DR/820, PC II
Merilno območje:	0 - 100 % Inhibition
Opis:	Toxicity Reagent set 25-49 tests
Parameter:	Toxicity
Številka metode:	10017
Število testov:	25 - 49
Vsebina paketa:	The ToxTrak Toxicity Reagent Set includes: Total Bacteria Count media set (15 TBC tubes, 15 individually wrapped sterile pipets), ToxTrak powder pillows (pk/50), ToxTrak accelerator solution (15mL), disposable culture tubes (5 packs; each pk/6), sterile transfer pipets (pk/15) and sodium thiosulfate standard (100 mL, 0.0246N).

Vsebina paketa

The ToxTrak Toxicity Reagent Set includes: Total Bacteria Count media set (15 TBC tubes, 15 individually wrapped sterile pipets), ToxTrak powder pillows (pk/50), ToxTrak accelerator solution (15mL), disposable culture tubes (5 packs; each pk/6), sterile transfer pipets (pk/15) and sodium thiosulfate standard (100 mL, 0.0246N).