



Laboratory Bottles, Polypropylene, Autoclavable, Wide Mouth, 500 mL, 12/pk

Številka izdelka: 1472459

EUR Cena: Kontakt
Odprema v 2 tednih

Polypropylene bottles suitable for sample collection to chemical storage. Matched polypropylene caps have a valve-seal design for a perfect, leakproof fit. A protruding ring under the cap holds shrink-wrapping securely in place. Autoclavable 121 °C (Do not autoclave with caps in place). Available in Narrow Mouth or Wide Mouth.

Specifikacije

Body Material:	Plastic
Kapaciteta:	500 mL
Material:	PP
Package size UOM:	12pk
Prostornina / velikost paketa:	12 /pk