



Test kit, acidity, model AC-DT

Številka izdelka:

2064000

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nearno



Acidity is usually caused by the presence of mineral acids, salts of strong acids, or free carbon dioxide. In surface or ground waters, acidity may form natural substances or industrial pollution. Excessive levels cause corrosion and can be detrimental to fish. Utilizing the Hach Digital Titrator, the Model AC-DT kit is used to determine acidity concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10-4,000 mg/L. Contains approximately 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Specifikacije

Ime metode:	Digital Titrator/Sodium hydroxide
Komentar:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Merilno območje:	10 - 4,000 mg/L CaCO ₃
Najmanjši korak rasti:	0.1-1.10
Parameter:	Acidity - (P or Total) as CaCO ₃
Platforma:	Digitalni Titrator
Število testov:	100
Test Kit model:	AC-DT
Vrsta primera:	F