



Be Right™



Test kit, chloride, model CD-51

Številka izdelka:

208601

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nevorno



Found in nearly all natural waters, chloride affects human taste above 250 mg/L. High levels inhibit plant growth, and many industrial processes need to limit chloride concentration. The model CD-DT, multi-range kit, is supplied with the demineralizer bottle for high-range sample dilution. Contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Ime metode:

Drop count titration/

Silver nitrate

Komentar:

*mg/L unless otherwise noted, ppb = $\mu\text{g}/\text{L}$, ppm = mg/L; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.

Merilno območje:

500 - 10000 mg/L & 5000 - 100000 mg/L Cl^-

Parameter:

Chloride - As Cl^-

Platforma:

Drop Count Titration

Smallest increment:

500

5000

Število testov:

100

Test Kit model:

CD-51

Vrsta primera:

B