



Cylinder, Graduated, Mixing, Glass, 10 mL ± 0.1 mL, 0.2 mL divisions, glass stopper #13

Številka izdelka: 2088638

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

This glass graduated cylinder is calibrated “to contain” with a molded in base. Each includes glass penny head stopper. Graduations are in durable white enamel bake.

Specifikacije

Certificiranje:	Class B
Kapaciteta:	10 mL
Material ohišja:	Glass
Merilna metoda:	To Contain
Prostornina / velikost paketa:	1 each
Vrsta:	eno metričen, mešan