



Be Right™



Centrifuge Tube, Conical Bottom, 10/pk

Številka izdelka: 2115839

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Centrifuge tube with conical bottom. Polycarbonate with molded-in 0.1mL graduations. For low-speed centrifugation up to 6000 x g. Autoclavable. The 15 mL size is also used with Settlometer Kit (P/N 21157-00). Supplied in packages of 10, without closures. Autoclavable. Nalge.

Specifikacije

Dolžina:	120 mm
Kapaciteta:	15 mL
Prostornina / velikost paketa:	10 /pk