



Test kit, chromium, model CH-14

Številka izdelka: 222702
EUR Cena: Kontakt
Datum pošiljanja ni prikazan

Nearno



Widely used as a corrosion inhibitor in cooling towers and engine cooling systems, chromium is also used in surface finishing, leather tanning, and other industries. The most widely used form, hexavalent chromium, is very toxic and strict limits are placed on discharges. Utilizing the drop count titration method to determine chromium (high range), this kit contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Ime metode:	Drop count titration/ Thiosulfate
Merilno območje:	5 - 100 mg/L & 50 - 1000 mg/L Cr ⁶⁺
Parameter:	Chromium, high range - As Hexavalent Cr
Platforma:	Drop Count Titration
Smallest increment:	5 50
Število testov:	100
Test Kit model:	CH-14
Vrsta primera:	D