



SA filter, cone, diameter 110 mm, soil analysis grade, 500/pk

Številka izdelka:	2263356
EUR Cena:	Kontakt
Datum pošiljanja ni prikazan	

Designed for sandy soils, SA filters are conically folded and have been optimized for use with the common extraction solutions used in soil nutrient analysis (Mehlich 1, Mehlich 3, Morgan, Modified Morgan, Bray, and Olsen). SA Filters provide the best combination of retention rating, filtration speed and purity requirements for soil nutrient testing (P, K, Ca, Mg). Medium filter retention.

Specifikacije

Diameter:	110 mm
Material:	Cellulose
Prostornina / velikost paketa:	50 /pk