



Test kit, hydrogenperoxide, model HYP-1

Številka izdelka: 2291700

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Nearno



Model HYP-1 Hydrogen Peroxide Test Kit (as H_2O_2) utilizes drop count titration and contains approximately 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Ime metode:	Drop count titration/ Thiosulfate
Merilno območje:	0,2 - 2,0 mg/L & 1 - 10 mg/L
Parameter:	Hydrogen Peroxide - As H_2O_2
Platforma:	Drop Count Titration
Smallest increment:	0.2
	1
Število testov:	100
Test Kit model:	HYP-1
Vrsta primera:	D