



Be Right™



Test kit, ascorbic acid, model ASC-1

Številka izdelka:

2308100

Nevarno

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan



Model ASC-1 Ascorbic Acid Test Kit was developed in cooperation with beverage producers and uses the drop-count titration method to measure ascorbic acid (vitamin C). Utilizes drop count titration with a range of 10-200 mg/L. Contains approximately 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Ime metode: Drop count titration/

Iodometric

Komentar: *mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.

Merilno območje: 10 - 200 mg/L

Parameter: Ascorbic Acid

Platforma: Drop Count Titration

Smallest increment: 10

Število testov: 100

Test Kit model: ASC-1

Vrsta primera: D