



Filter, membrane, black, pore size 0.45 µm, diameter 47 mm, 200/pk, w/out pads

Številka izdelka: 2490300

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Provides contrast for microbiological and particulate analysis. Modified polyethersulfone creates dark background for excellent contrast useful in enumerating opaque colonies and light-covered particulate. White grid line aids enumeration. Certified for membrane filtration analysis. Sterile and individually packaged. Pore size 0.45 µm. Diameter: 47 mm.

Specifikacije

Diameter:	47 mm
Filter Type:	Filter
Material:	Modified polyethersulfone
Prostornina / velikost paketa:	200 /pk
Velikost por:	0.45 µm