



Test tube/ MPN rack, polyethylene, 20 mm (O.D.), 40 holes

Številka izdelka: 2497902

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Forty-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5").

Specifikacije

Asortiman testov:	4 x 10
Diameter:	20 mm
Luknje za teste:	40
Prostornina / velikost paketa:	1 each