



Test tube, rack, polyethylene, 25 mm (O.D.), 24 holes

Številka izdelka: 2497903

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Twenty to ninety-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5").

Specifikacije

Asortiman testov:	3 x 8
Diameter:	25 mm
Luknje za teste:	24
Prostornina / velikost paketa:	1 each