



Test kit, chloride, model CDS-DT

Številka izdelka: 2580600
EUR Cena: Kontakt
Datum pošiljanja ni prikazan

Nearno



The single-parameter Model CDS-DT uses the Hach Digital Titrator with Silver Nitrate Titrant Cartridges to determine chloride concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10-10,000 mg/L. Contains approximately 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Specifikacije

Ime metode:	Digital Titrator/ Silver nitrate
Merilno območje:	10 - 10,000 mg/L Cl ⁻
Parameter:	Chloride - As Cl ⁻ , 10-10,000 mg/L Cl ⁻
Platforma:	Digitalni Titrator
Smallest increment:	10
Število testov:	100
Test Kit model:	CDS-DT
Vrsta primera:	F

Obvezni pripomočki

- Pipeta, spremenljiva prostornina; 1,0 - 5,0 mL (Item BBP065)