



Test kit, hypochlorite, model CN-HR

Številka izdelka:

2687200

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nearno



Hypochlorite Kit, Model CN-HR - Drop Count

Utilizing the simple, economical drop count titration method, this kit is used to determine high range hypochlorite solution concentrations (bleach solutions). Range: 5-15%, contains 100 tests.

Simple, economical titration method

Kit contains all required reagents, solutions, and glassware

Specifikacije

Ime metode:	Drop count titration/ Thiosulfate
Merilno območje:	50 - 150 g/L Cl_2 (5-15%)
Parameter:	Hypochlorite, high range - (Bleach) As Cl_2
Platforma:	Drop Count Titration
Smallest increment:	0.5%
Število testov:	100
Test Kit model:	CN-HR
Vrsta primera:	D