



## APA6000 Alkalinity Standard 2, 500 mg/L, 1 L

Številka izdelka:	2826253
EUR Cena:	Kontakt
Datum pošiljanja ni prikazan	

The APA6000 Alkalinity Process analyzer uses five reagents each of which is supplied in a 1 L bottle designed to provide reagents for 30 days of continuous use. APA Alkalinity Reagent 1 (2826153) is an acid titrant that reacts with the alkalinity in the sample. APA Alkalinity Reagent 2 (2696653) is an indicator used to determine the endpoints of the titration. APA Alkalinity Standards 1 and 2 (6001100) are used together to perform a multiple point calibration. The fifth reagent, APA Alkalinity Cleaning Solution (2697053), cleans the system during the Prime and Instrument Clean functions.

---

### Specifikacije

Inštrument:	APA6000, EZ4004
Koncentracija:	500 mg/L
Parameter:	Alkalinity
Pogoji skladiščenja:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.