



Filling Solution, Reference, 3 M KCl with AgCl, 28 mL

Številka izdelka: 2841700

EUR Cena: Kontakt

Na voljo

Filling solution for electrode reference element, Potassium Chloride (KCl) with Saturated AgCl, 3 M, 28 mL bottle

Specifikacije

Formula:	KCl 3M & AgCl Sat.
Opis:	MTC301, PHC301, 5012T
Platforma:	Raztopina
Prostornina / velikost paketa:	28 mL
Solution Type:	Refill